

CLEAR *Light*[™] and **Clear100**[™] Acne PhotoClearing[™] Systems

Rapid Results with No Side Effects

ClearLight and Clear100 Selective Photo Clearing (SPC) systems have a clearance rate of 60%–70% of inflammatory lesions in up to one-third the time of traditional methods and without side effects.

Use alone or in conjunction with topical therapies to offer long-lasting clearance.



- **Innovative new technology**

Unique light-based Selective Photo Clearing (SPC) technology precisely destroys surface and ductal p. acnes.

- **Faster treatment**

Up to three times faster than standard topical or oral treatments. Early visible results encourage patient compliance.

- **Hands-free**

Non-contact treatment saves time.

- **No side effects**


- **High patient throughput**

Rapid 20-minute, full-face treatments can be administered by staff.




CLEAR*Light*™ and Clear100™

Specifications		ClearLight	Clear100
Technology	SPC System, double head	X	X
Wavelength	405-420 nm	X	X
Energy Density	60 Joules/cm ²	X	X
Treatment Area	2 x 30 x 30 cm	X	X
Control and Display	Computerized, clear LCD display	X	X
Treatment Ergonomic	One-touch, hands-free treatment; electrical one-touch treatment head positioning	X	X
Electrical Requirements	220V; 16A; 50 Hz 110V; 20A; 60 Hz	X	X
Size	50 x 41 x 162 cm / 21 x 16 x 64 in (WDH) (folded) 131 x 40 x 163 cm / 52 x 16 x 64 in (WDH) (open)	X	X
Weight	80 kg / 176 lbs	X	X
Integrated Patient Database		X	
Integrated Digital Camera		X	



Acne Selective Photo Clearing (SPC) technology vs. other treatment modalities as based on a multi-center study in USA, Spain, Israel and Belgium. Preliminary studies show a non-response rate of approximately 15%.




Selective light safely and accurately penetrates the skin and destroys acne bacteria without causing harm to tissue.


Before and After



Photos: Hector Leal-Silva, MD



Photos: Yoram Harth, MD



Lumenis Global Sales Offices

Lumenis Inc.
2400 Condensa Street
Santa Clara, CA 95051 USA
Tel 408 764 3000
877 586 3647
Fax 408 764 3999
800 505 1133

Lumenis (HOLLAND) BV
Amstelveen, The Netherlands
Tel +31 20 347 50 60
Fax +31 20 347 50 99

Lumenis (GERMANY) GmbH
Dreieich-Dreieichenhain,
Germany
Tel +49 6071 922 0
Fax +49 6071 922 400

Lumenis (FRANCE) SARL
Orsay, France
Tel +33 1 69 33 14 20
Fax +33 1 60 19 57 23

Lumenis (ITALY) Srl.
Rome, Italy
Tel +39 06 90 75 230
Fax +39 06 90 75 269

Lumenis (UK) Ltd.
London, UK
Tel +44 20 8324 4200
Fax +44 20 8324 4222


Lumenis (JAPAN) Co. Ltd.
Tokyo, Japan
Tel +81 3 5789 8300
Fax +81 3 5789 8310

Lumenis (Holland) BV
Beijing Office, Beijing, China
Tel +86 10 6510 2620
Fax +86 10 6510 2621

Lumenis (HONG KONG) Ltd.
Hong Kong
Tel +852 2174 2800
Fax +852 2722 5151

Internet: www.lumenis.com

Email: information@lumenis.com



Lumenis, its logo, ClearLight, Clear100 and SPC are trademarks or registered trademarks of the Lumenis Group of Companies. Product specifications are subject to change without notice. ©2004, the Lumenis Group of Companies. All rights reserved.

Printed in the U.S.A.

Web only



A-225 10/04 REGISTERED